

Anti-FcRH5 Antibody (DFRF4539A Naked Antibody)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-FcRH5 Antibody (MFRF3266A, DFRF4539A Naked Antibody) is a highly specific monoclonal antibody targeting FcRH5 (also known as IRTA2 or CD307e) and serves as the antibody backbone for the ADC DFRF4539A (RG-7598), which is used in research related to relapsed or refractory multiple myeloma.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481