

Anti-TPBG Antibody (PF-06263507)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-TPBG Antibody (PF-06263507) is a human antibody expressed in CHO cells, targeting TPBG. It features a hulgG1 type heavy chain and a huk type light chain, with a predicted molecular weight (MW) of 145.8 kDa. For isotype control, refer to HumanIgG1kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481