


## Ginsenoside Rd2

## Chemical Properties

CAS No. :	83480-64-2
Formula:	C <sub>47</sub> H <sub>80</sub> O <sub>17</sub>
Molecular Weight:	917.14
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ginsenoside Rd2 (Quinquenoside L10) is a natural product with Nourish the strong, calm the mind, nourish the body.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0903 mL	5.4517 mL	10.9035 mL
5 mM	0.2181 mL	1.0903 mL	2.1807 mL
10 mM	0.109 mL	0.5452 mL	1.0903 mL
50 mM	0.0218 mL	0.109 mL	0.2181 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481