

## Compound T128647

## Chemical Properties

CAS No. :

Formula: C<sub>36</sub>H<sub>58</sub>O<sub>10</sub>

Molecular Weight: 650.85

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	3-{3-[(4E)-2,6-dihydroxy-6-methylhept-4-en-2-yl]-6,9a,9b-trimethyl-7-(3-oxoprop-1-en-2-yl)-4-[[[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy]-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl}propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5365 mL	7.6823 mL	15.3645 mL
5 mM	0.3073 mL	1.5365 mL	3.0729 mL
10 mM	0.1536 mL	0.7682 mL	1.5365 mL
50 mM	0.0307 mL	0.1536 mL	0.3073 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481