

NR 7h

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	NR 7h is a potent and selective degrader of p38 α and p38 β , exhibiting DC50 values of 24 nM for p38 α in T47D cells and 27.2 nM in MB-MDA-231 cells, respectively [1].
Targets(IC50)	p38 MAPK,PROTACs

Solubility Information

Solubility	DMSO: 99.89 mg/mL,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481