

## Anti-CD207 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

|               |   |
|---------------|---|
| Description   | Anti-CD207 antibody is a kind of antibody that can be used for activity assays, functional studies, binding affinity analysis, and other applications, serving as an essential tool in biological research. |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481