

## Rat IgG2a Isotype Control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Rat IgG2a isotype control is a rat IgG2a isotype-matched control antibody directed against trinitrophenol, and it is specifically designed for use as a negative control in immunological assays, including flow cytometry and immunohistochemistry, to validate antibody specificity and ensure accurate interpretation of experimental immune-response data.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481