

Anti-TNFR2/CD120b/TNFR1B Antibody (6I44)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-TNFR2/CD120b/TNFR1B Antibody (6I44) is a Rabbit IgG monoclonal antibody targeting Mouse TNFR2/CD120b/TNFR1B.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481