

## Anti-IL-23 P19/IL23A Antibody (8P117)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-IL-23 P19/IL23A Antibody (8P117) is a Mouse IgG1 monoclonal antibody targeting Human, Mouse IL-23 P19/IL23A.
Targets(IC50)	IL Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481