

## Anti-IL-23 P19/IL23A Antibody (7N865)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-IL-23 P19/IL23A Antibody (7N865) is a mouse (variable region) / human (kappa / IgG1 constant) chimeric monoclonal antibody targeting Human, Mouse IL-23 P19/IL23A.
Targets(IC50)	IL Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481