

## Anti-Enterovirus 71 VP4 Antibody (3E202)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Enterovirus 71 VP4 Antibody (3E202) is a monoclonal antibody targeting EV71 (Enterovirus 71) VP4, which is a mouse-derived antibody.
Targets(IC50)	Virus Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481