

Anti-SARS-CoV/SARS-CoV-2 Spike Antibody (2J777)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Anti-SARS-CoV/SARS-CoV-2 Spike Antibody (2J777) is a mouse (variable region) / human (kappa / IgG1 constant) chimeric monoclonal antibody targeting SARS SARS-CoV/SARS-CoV-2 Spike. |
| Targets(IC50) | SARS-CoV |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481