

Anti-Influenza A H10 Hemagglutinin/HA Antibody (8B805)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Influenza A H10 Hemagglutinin/HA Antibody (8B805) is a Mouse IgG1 monoclonal antibody targeting Influenza Influenza A H10 Hemagglutinin/HA.
Targets(IC50)	Influenza Virus

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481