

## Anti-SARS-CoV-2 Omicron (BA.4/BA.5/BA.2.75) Spike RBD Specific Antibody (3M568)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-SARS-CoV-2 Omicron (BA.4/BA.5/BA.2.75) Spike RBD Specific Antibody (3M568) is a mouse (variable region) / human (kappa / IgG1 constant) chimeric monoclonal antibody targeting SARS-CoV-2 SARS-CoV-2 Omicron (BA.4/BA.5/BA.2.75) Spike RBD
Targets(IC50)	SARS-CoV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481