

Aseanostatin P5

Chemical Properties


CAS No. : 5502-94-3

Formula: C₁₅H₃₀O₂

Molecular Weight: 242.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Aseanostatin P5 inhibits myeloperoxidase (MPO) release from human polymorphonuclear leukocytes.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1254 mL	20.6271 mL	41.2541 mL
5 mM	0.8251 mL	4.1254 mL	8.2508 mL
10 mM	0.4125 mL	2.0627 mL	4.1254 mL
50 mM	0.0825 mL	0.4125 mL	0.8251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481