


Thalifoline

Chemical Properties

CAS No. :	21796-15-6
Formula:	C ₁₁ H ₁₃ NO ₃
Molecular Weight:	207.229
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Thalifoline shows antifungal activity, with 100 ug against Fusarium oxysporum f. sp. lycopersici.
Targets(IC50)	Antifection,Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8256 mL	24.1278 mL	48.2556 mL
5 mM	0.9651 mL	4.8256 mL	9.6511 mL
10 mM	0.4826 mL	2.4128 mL	4.8256 mL
50 mM	0.0965 mL	0.4826 mL	0.9651 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481