

## Compound T129042

### Chemical Properties


CAS No. :

Formula: C42H68O13

Molecular Weight: 780.99

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	7-[[[(2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-[[[(2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy]oxan-2-yl]oxy]-1-[(4E)-2,6-dihydroxy-6-methylhept-4-en-2-yl]-3a,6,6,9b,11a-pentamethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-10-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2804 mL	6.4021 mL	12.8043 mL
5 mM	0.2561 mL	1.2804 mL	2.5609 mL
10 mM	0.128 mL	0.6402 mL	1.2804 mL
50 mM	0.0256 mL	0.128 mL	0.2561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481