

## 2-(hydroxymethyl)-6-[(2-methylbut-3-en-2-yl)oxy]oxane-3,4,5-triol

## Chemical Properties

CAS No. :

Formula: C<sub>11</sub>H<sub>20</sub>O<sub>6</sub>

Molecular Weight: 248.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	2-(hydroxymethyl)-6-[(2-methylbut-3-en-2-yl)oxy]oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0277 mL	20.1386 mL	40.2771 mL
5 mM	0.8055 mL	4.0277 mL	8.0554 mL
10 mM	0.4028 mL	2.0139 mL	4.0277 mL
50 mM	0.0806 mL	0.4028 mL	0.8055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481