


Ganoderic acid E

Chemical Properties

CAS No. : 98665-14-6
 Formula: C₃₀H₄₀O₇
 Molecular Weight: 512.63
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ganoderic acid F has anti-hepatitis B, anti-inflammatory, and anti-tumor-promoting activities.
Targets(IC50)	Antiviral

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9507 mL	9.7536 mL	19.5072 mL
5 mM	0.3901 mL	1.9507 mL	3.9014 mL
10 mM	0.1951 mL	0.9754 mL	1.9507 mL
50 mM	0.039 mL	0.1951 mL	0.3901 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481