


AS1134900

## Chemical Properties

CAS No. : 149760-98-5  
 Formula: C<sub>25</sub>H<sub>22</sub>N<sub>8</sub>S  
 Molecular Weight: 466.56  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	AS1134900 is a highly selective, non-competitive, metabotropic NADP <sup>+</sup> -dependent malic enzyme 1 (ME1) inhibitor.
Targets(IC50)	Others,Wnt/beta-catenin

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1433 mL	10.7167 mL	21.4335 mL
5 mM	0.4287 mL	2.1433 mL	4.2867 mL
10 mM	0.2143 mL	1.0717 mL	2.1433 mL
50 mM	0.0429 mL	0.2143 mL	0.4287 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481