

Compound T128652

Chemical Properties

CAS No. :

Formula: C₄₁H₇₀O₁₂

Molecular Weight: 755

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2R,3R,4S,5S,6R)-2-({1-[(4E)-2,6-dihydroxy-6-methylhept-4-en-2-yl]-3a,3b,6,6,9a-pentamethyl-7-[(2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy}-hexadecahydro-1H-cyclopenta[a]phenanthren-11-yl}oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3245 mL	6.6225 mL	13.245 mL
5 mM	0.2649 mL	1.3245 mL	2.649 mL
10 mM	0.1325 mL	0.6623 mL	1.3245 mL
50 mM	0.0265 mL	0.1325 mL	0.2649 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481