

## Compound T130503

### Chemical Properties


CAS No. :

Formula: C<sub>25</sub>H<sub>26</sub>O<sub>8</sub>

Molecular Weight: 454.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	methyl 3-(3,4-dihydroxyphenyl)-2-[[4-hydroxy-3-(3-methylbut-2-en-1-yl)phenyl]methyl]-4-methoxy-5-oxo-2,5-dihydrofuran-2-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2003 mL	11.0016 mL	22.0032 mL
5 mM	0.4401 mL	2.2003 mL	4.4006 mL
10 mM	0.220 mL	1.1002 mL	2.2003 mL
50 mM	0.044 mL	0.220 mL	0.4401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481