

2-[5-(2-hydroxypropyl)oxolan-2-yl]propanoic acid

Chemical Properties

CAS No. :

Formula: C₁₀H₁₈O₄

Molecular Weight: 202.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-[5-(2-hydroxypropyl)oxolan-2-yl]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9444 mL	24.7219 mL	49.4438 mL
5 mM	0.9889 mL	4.9444 mL	9.8888 mL
10 mM	0.4944 mL	2.4722 mL	4.9444 mL
50 mM	0.0989 mL	0.4944 mL	0.9889 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481