


## Goyaglycoside d

## Chemical Properties

CAS No. : 333332-50-6  
 Formula: C38H62O9  
 Molecular Weight: 662.905  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Goyaglycoside d is a natural product for research related to life sciences. The catalog number is TN6247 and the CAS number is 333332-50-6.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5085 mL	7.5425 mL	15.085 mL
5 mM	0.3017 mL	1.5085 mL	3.017 mL
10 mM	0.1509 mL	0.7543 mL	1.5085 mL
50 mM	0.0302 mL	0.1509 mL	0.3017 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481