

Compound 0080-0033

## Chemical Properties


CAS No. :

Formula: C<sub>27</sub>H<sub>38</sub>N<sub>2</sub>O<sub>7</sub>

Molecular Weight: 502.608

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0080-0033 is a useful organic compound for research related to life sciences and the catalog number is T131699.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9896 mL	9.9481 mL	19.8961 mL
5 mM	0.3979 mL	1.9896 mL	3.9792 mL
10 mM	0.199 mL	0.9948 mL	1.9896 mL
50 mM	0.0398 mL	0.199 mL	0.3979 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481