


## 5-Iodocytosine

## Chemical Properties

CAS No. :	1122-44-7
Formula:	C <sub>4</sub> H <sub>4</sub> IN <sub>3</sub> O
Molecular Weight:	237
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-Iodocytosine is a Heterocyclic compound-pyrimidine, intermediate and building block, halo-.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2194 mL	21.097 mL	42.1941 mL
5 mM	0.8439 mL	4.2194 mL	8.4388 mL
10 mM	0.4219 mL	2.1097 mL	4.2194 mL
50 mM	0.0844 mL	0.4219 mL	0.8439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481