

3,3'-[Iminobis(methylene)]bis-2(3H)furanone

Chemical Properties


CAS No. :

Formula: C₁₀H₁₅NO₄

Molecular Weight: 213.233

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,3'-[Iminobis(methylene)]bis-2(3H)furanone is a useful organic compound for research related to life sciences and the catalog number is T123845.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6898 mL	23.4489 mL	46.8977 mL
5 mM	0.938 mL	4.6898 mL	9.3795 mL
10 mM	0.469 mL	2.3449 mL	4.6898 mL
50 mM	0.0938 mL	0.469 mL	0.938 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481