


Pipemidic acid

Chemical Properties

CAS No. :	51940-44-4
Formula:	C ₁₄ H ₁₇ N ₅ O ₃
Molecular Weight:	303.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pipemidic acid (Acido pipemidico), a new antibacterial agent, inhibits the activity of Pseudomonas aeruginosa.
Targets(IC50)	Antibacterial, Antibiotic

Solubility Information

Solubility	DMSO: < 1 mg/mL (insoluble or slightly soluble), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2968 mL	16.4842 mL	32.9685 mL
5 mM	0.6594 mL	3.2968 mL	6.5937 mL
10 mM	0.3297 mL	1.6484 mL	3.2968 mL
50 mM	0.0659 mL	0.3297 mL	0.6594 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Hoeglund AB, et al. J Med Chem. 2010 Feb 11;53(3):1056-66.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481