


Beta-D-Glucopyranosyl nitromethane

Chemical Properties

CAS No. :	81846-60-8
Formula:	C7H13NO7
Molecular Weight:	223.18
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Beta-D-Glucopyranosyl nitromethane is a glycosyl derivative of nitromethane, functioning as a salt of a highly basic anion exchanger in the OH cycle.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4807 mL	22.4034 mL	44.8069 mL
5 mM	0.8961 mL	4.4807 mL	8.9614 mL
10 mM	0.4481 mL	2.2403 mL	4.4807 mL
50 mM	0.0896 mL	0.4481 mL	0.8961 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Petrus, L., et al. Preparation of some glycosyl derivatives of nitromethane. *Chemicke Zvesti* (1982), 36(1), 103-110.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481