

Compound T128714

Chemical Properties

CAS No. :

Formula: C₂₆H₃₆O₇

Molecular Weight: 460.57

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	methyl 2,5-dihydroxy-2,3b,6,6,9a,11a-hexamethyl-11-methylidene-1,3,7-trioxo-hexadecahydro-1H-cyclopenta[a]phenanthrene-3a-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1712 mL	10.8561 mL	21.7122 mL
5 mM	0.4342 mL	2.1712 mL	4.3424 mL
10 mM	0.2171 mL	1.0856 mL	2.1712 mL
50 mM	0.0434 mL	0.2171 mL	0.4342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481