

EIDD-1619

Chemical Properties


CAS No. : 1659302-92-7

Formula: C₂₂H₃₄NNa₂O₆P

Molecular Weight: 485.468

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	EIDD-1619 is a novel water-soluble analog of progesterone.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0599 mL	10.2993 mL	20.5986 mL
5 mM	0.412 mL	2.0599 mL	4.1197 mL
10 mM	0.206 mL	1.0299 mL	2.0599 mL
50 mM	0.0412 mL	0.206 mL	0.412 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481