

Boc-N-PEG1-C2-NHS ester

Chemical Properties

CAS No. : 1260092-55-4


Formula: C₁₄H₂₂N₂O₇

Molecular Weight: 330.33

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boc-N-PEG1-C2-NHS ester, a PEG-based linker, finds utility in PROTAC synthesis[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands linked together: one targets an E3 ubiquitin ligase, and the other binds to the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0273 mL	15.1364 mL	30.2728 mL
5 mM	0.6055 mL	3.0273 mL	6.0546 mL
10 mM	0.3027 mL	1.5136 mL	3.0273 mL
50 mM	0.0605 mL	0.3027 mL	0.6055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481