


LY-393615

Chemical Properties

CAS No. : 325819-97-4
 Formula: C₂₁H₂₆ClF₂NO
 Molecular Weight: 381.89
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	LY-393615 is a neuronal calcium channel blocker.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6186 mL	13.0928 mL	26.1856 mL
5 mM	0.5237 mL	2.6186 mL	5.2371 mL
10 mM	0.2619 mL	1.3093 mL	2.6186 mL
50 mM	0.0524 mL	0.2619 mL	0.5237 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481