

MR837

Chemical Properties


CAS No. : 1210906-48-1

Formula: C₁₆H₁₄N₂O₅

Molecular Weight: 282.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	MR837 (NSD2-PWWP1 antagonist 3f) is a NSD2-PWWP1 antagonist.
Targets(IC50)	Histone Methyltransferase

Solubility Information

Solubility	DMSO: 11 mg/mL (38.96 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (7.08 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5416 mL	17.7079 mL	35.4158 mL
5 mM	0.7083 mL	3.5416 mL	7.0832 mL
10 mM	0.3542 mL	1.7708 mL	3.5416 mL
50 mM	0.0708 mL	0.3542 mL	0.7083 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481