

## m-PEG20-alcohol

## Chemical Properties

CAS No. :	1059605-06-9
Formula:	C41H84O21
Molecular Weight:	913.102
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	m-PEG20-alcohol is a PEG-based linker for PROTACs, joining two essential ligands crucial for forming PROTAC molecules and enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are bifunctional molecules composed of two ligands linked together: one binding to an E3 ubiquitin ligase and the other targeting a specific protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0952 mL	5.4759 mL	10.9517 mL
5 mM	0.219 mL	1.0952 mL	2.1903 mL
10 mM	0.1095 mL	0.5476 mL	1.0952 mL
50 mM	0.0219 mL	0.1095 mL	0.219 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481