

Compound N105-0100

Chemical Properties


CAS No. :

Formula: C₂₆H₄₄O₄

Molecular Weight: 420.634

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N105-0100 is a useful organic compound for research related to life sciences and the catalog number is T131146.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3774 mL	11.8869 mL	23.7739 mL
5 mM	0.4755 mL	2.3774 mL	4.7548 mL
10 mM	0.2377 mL	1.1887 mL	2.3774 mL
50 mM	0.0475 mL	0.2377 mL	0.4755 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481