


## Brilliant Blue R250

## Chemical Properties

CAS No. :	6104-59-2
Formula:	C45H45N3NaO7S2
Molecular Weight:	826.97
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	CBBR is an inhibitor of wild type and mutant alpha-synuclein aggregation. CBBR also a modulator of neurotoxicity.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: Soluble, ( $< 1$ mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2092 mL	6.0462 mL	12.0923 mL
5 mM	0.2418 mL	1.2092 mL	2.4185 mL
10 mM	0.1209 mL	0.6046 mL	1.2092 mL
50 mM	0.0242 mL	0.1209 mL	0.2418 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

### Reference

- Ahsan N, Siddique IA, Gupta S, Surolia A. A routinely used protein staining dye acts as an inhibitor of wild type and mutant alpha-synuclein aggregation and modulator of neurotoxicity. *Eur J Med Chem.* 2018 Jan 1;143:1174-1184. doi: 10.1016/j.ejmech.2017.10.002. Epub 2017 Oct 12. PubMed PMID: 29150334.
- Yari S, Hadizadeh Tasbiti A, Ghanei M, Siadat SD, Yari F, Bahrmann A. Proteome-scale MDR-TB-antibody responses for identification of putative biomarkers for the diagnosis of drug-resistant Mycobacterium tuberculosis. *Int J Mycobacteriol.* 2016 Dec;5 Suppl 1:S134-S135. doi: 10.1016/j.ijmyco.2016.10.011. Epub 2016 Oct 27. PubMed PMID: 28043508.
- Jaiswal AK, Sudan V, Pandey V, Singh A, Gaur RS, Kanojiya D, Nigam R, Shanker D. Sex dependent alterations in the protein characterization patterns of *Haemonchus contortus*. *J Parasit Dis.* 2016 Sep;40(3):1006-8. doi: 10.1007/s12639-014-0624-2. Epub 2014 Dec 20. PubMed PMID: 27605828; PubMed Central PMCID: PMC4996237.
- Yasumitsu H, Ozeki Y, Kanaly RA. RAMA casein zymography: Time-saving and highly sensitive casein zymography for MMP7 and trypsin. *Electrophoresis.* 2016 Nov;37(22):2959-2962. doi: 10.1002/elps.201600346. Epub 2016 Sep 9. PubMed PMID: 27500714.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481