


D&C Red No. 8

Chemical Properties

CAS No. :	2092-56-0
Formula:	C ₁₇ H ₁₃ ClN ₂ NaO ₄ S
Molecular Weight:	399.8
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	D&C Red No. 8 is a dye.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5013 mL	12.5063 mL	25.0125 mL
5 mM	0.5003 mL	2.5013 mL	5.0025 mL
10 mM	0.2501 mL	1.2506 mL	2.5013 mL
50 mM	0.050 mL	0.2501 mL	0.5003 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481