


(S)-Ketorolac

Chemical Properties

CAS No. :	66635-92-5
Formula:	C ₁₅ H ₁₃ NO ₃
Molecular Weight:	255.27
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-Ketorolac is a potent inhibitor of COX1 and COX2 enzyme, and is nonsteroidal anti-inflammatory agent.
Targets(IC50)	Others,COX

Solubility Information

Solubility	DMSO: 200 mg/mL (783.48 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9174 mL	19.5871 mL	39.1742 mL
5 mM	0.7835 mL	3.9174 mL	7.8348 mL
10 mM	0.3917 mL	1.9587 mL	3.9174 mL
50 mM	0.0783 mL	0.3917 mL	0.7835 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Handley DA, et al. Preclinical enantioselective pharmacology of (R)- and (S)-ketorolac. J Clin Pharmacol. 1998 Feb; 38(2S):25S-35S.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481