

K 309

## Chemical Properties


CAS No. : 52321-54-7

Formula: C<sub>16</sub>H<sub>13</sub>NO<sub>2</sub>S

Molecular Weight: 283.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	K 309 is a bioactive chemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5292 mL	17.646 mL	35.292 mL
5 mM	0.7058 mL	3.5292 mL	7.0584 mL
10 mM	0.3529 mL	1.7646 mL	3.5292 mL
50 mM	0.0706 mL	0.3529 mL	0.7058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481