


Lipoamido-PEG4-acid

Chemical Properties

CAS No. : 1314378-10-3

Formula: C19H35NO7S2

Molecular Weight: 453.61



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lipoamido-PEG4-acid is a PEG-based linker for PROTACs that joins two essential ligands, crucial for forming PROTAC molecules, enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one binds to an E3 ubiquitin ligase, and the other targets a specific protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2045 mL	11.0227 mL	22.0454 mL
5 mM	0.4409 mL	2.2045 mL	4.4091 mL
10 mM	0.2205 mL	1.1023 mL	2.2045 mL
50 mM	0.0441 mL	0.2205 mL	0.4409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481