

## Cymen-ar-ol

## Chemical Properties

CAS No. : 39660-61-2

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Cymen-ar-ol

## Biological Description

Description	Cymen-ar-ol is a broad spectrum and highly effective bactericide that can contact mucous membrane.
Targets(IC50)	Antibacterial

## Solubility Information

Solubility	DMSO: 60 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---------------------------------------------------------------------------------------------------------------

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481