


SDZ 21009

Chemical Properties

CAS No. : 39731-05-0
 Formula: C₁₉H₂₈N₂O₄
 Molecular Weight: 348.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	β -adrenoceptor and 5-HT _{1A/1B} receptor antagonist
Targets(IC ₅₀)	Others,5-HT Receptor

Solubility Information

Solubility	DMSO: <34.84 mg/mL,Sonication is recommended. Ethanol: <34.84 mg/mL,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8699 mL	14.3497 mL	28.6993 mL
5 mM	0.574 mL	2.8699 mL	5.7399 mL
10 mM	0.287 mL	1.435 mL	2.8699 mL
50 mM	0.0574 mL	0.287 mL	0.574 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481