


m-PEG11-Tos

Chemical Properties

CAS No. :	58320-73-3
Formula:	C28H50O13S
Molecular Weight:	
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	M-PEG11-Tos is a PEG derivative containing a tosyl group. The hydrophilic PEG spacer increases solubility in aqueous media. For nucleophilic substitution reactions, the tosyl group is a good leaving group.
Targets(IC50)	Others,PROTAC Linker

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Reference

- Brown A, Patel S, Ward C, Lorenz A, Ortiz M, DuRoss A, Wieghardt F, Esch A, Otten EG, Heiser LM, Korolchuk VI, Sun C, Sarkar S, Sahay G. PEG-lipid micelles enable cholesterol efflux in Niemann-Pick Type C1 disease-based lysosomal storage disorder. *Sci Rep.* 2016 Aug 30;6:31750. doi: 10.1038/srep31750. PubMed PMID: 27572704; PubMed Central PMCID: PMC5004151.
- Saifer MG, Williams LD, Sobczyk MA, Michaels SJ, Sherman MR. Selectivity of binding of PEGs and PEG-like oligomers to anti-PEG antibodies induced by methoxyPEG-proteins. *Mol Immunol.* 2014 Feb;57(2):236-46. doi: 10.1016/j.molimm.2013.07.014. Epub 2013 Nov 5. PubMed PMID: 24200843.
- Garay RP, El-Gewely R, Armstrong JK, Garratty G, Richette P. Antibodies against polyethylene glycol in healthy subjects and in patients treated with PEG-conjugated agents. *Expert Opin Drug Deliv.* 2012 Nov;9(11):1319-23. doi: 10.1517/17425247.2012.720969. Epub 2012 Aug 30. PubMed PMID: 22931049.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481