


Ivaltinostat

Chemical Properties

CAS No. :	936221-33-9
Formula:	C ₂₄ H ₃₃ N ₃ O ₄
Molecular Weight:	427.54
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	CG-200745 is a potent inhibitor of HDAC. CG200745 induces apoptosis and also inhibits the deacetylation of histone H3 and tubulin.
Targets(IC50)	Apoptosis,HDAC,MDM-2/p53
In vitro	CG200745 induces the accumulation of p53. It promotes p53-dependent transactivation, and it also enhances the expression of MDM2 and p21 (Waf1/Cip1) proteins [1,2]. CG-200745 (0-10 μM; 48 hours) reduces the Calu6 cells proliferation to 40% of untreated cells. CG-200745 (0-10 μM; 1-24 hours) treatment with low concentration significantly increases the acetylation of histone H3 and H4 in Calu6 cells at various sites in a time-dependent manner up to 24 hours after treatment. CG-200745 (3 μM; 1-24 hours) significantly increases Calu6 cells proportion in the G2/M phase (69%) [3].
In vivo	CG-200745 (p.o.; 30?mg/kg/day; for 7 days) weaken inflammatory cytokines, adhesion molecules, oxidative stress in UUO kidneys [4].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.339 mL	11.6948 mL	23.3896 mL
5 mM	0.4678 mL	2.339 mL	4.6779 mL
10 mM	0.2339 mL	1.1695 mL	2.339 mL
50 mM	0.0468 mL	0.2339 mL	0.4678 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Oh ET, et al. Novel histone deacetylase inhibitor CG200745 induces clonogenic cell death by modulating acetylation of p53 in cancer cells. Invest New Drugs. 2012 Apr;30(2):435-42.

Hwang JJ, et al. A novel histone deacetylase inhibitor, CG200745, potentiates anticancer effect of docetaxel in prostate cancer via decreasing Mcl-1 and Bcl-XL. Invest New Drugs. 2012 Aug;30(4):1434-42.

Chun SM, et al. Epigenetic modulation with HDAC inhibitor CG200745 induces anti-proliferation in non-small cell lung cancer cells. PLoS One. 2015 Mar 17;10(3):e0119379.

Choi HS, et al. Histone deacetylase inhibitor, CG200745 attenuates renal fibrosis in obstructive kidney disease. Sci Rep. 2018 Aug 1;8(1):11546.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481