


RepSox

Chemical Properties

CAS No. :	446859-33-2
Formula:	C17H13N5
Molecular Weight:	287.32
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	RepSox (ALK5 Inhibitor II) is a TGFβR-1/ALK5 inhibitor that selectively inhibits the binding of ATP to ALK5 and the autophosphorylation of ALK5 (IC50=23/4 nM). RepSox induces adipogenesis in MEFs cells.
Targets(IC50)	ALK,TGF-beta/Smad
In vitro	<p>METHODS: Mouse embryonic fibroblast MEFs were treated with RepSox (10 μM) for 10 days and cellular lipids were detected using oil red O staining.</p> <p>RESULTS: Many adipocyte-like cells appeared in RepSox-treated MEFs, and oil red O staining of lipid droplets confirmed adipocyte differentiation. [1]</p> <p>METHODS: Sheep fibroblasts were treated with RepSox (10-25 nM) for 3 days, and cell morphology was observed using a microscope.</p> <p>RESULTS: After three days of RepSox treatment, sheep fibroblasts acquired a significantly different morphology and changed from a spindle shape to an elongated shape with more clusters. [2]</p>
In vivo	<p>METHODS: To test the effect on the differentiation of EGCs from intestinal glial cells, RepSox (3-10 mg/kg in 0.5% CMC-Na) was injected intraperitoneally to GFAP-Cre: Rosa26-tdTomato mice once daily for two weeks.</p> <p>RESULTS: RepSox promotes the conversion of EGCs into neurons and improves gastrointestinal motility in vivo. [3]</p> <p>METHODS: To detect anti-tumor activity in vivo, RepSox (5-20 mg/kg) was injected intraperitoneally into a nude mouse model bearing osteosarcoma 143B every two days for three weeks.</p> <p>RESULTS: RepSox effectively inhibited the growth of osteosarcoma and had low toxicity in nude mice. [4]</p>
Kinase Assay	The kinase domain of ALK5 is cloned by PCR and expressed in a baculovirus/Sf9 cells system. The protein is 6-His tagged in the C terminus and purified by affinity chromatography using a Ni ²⁺ -column, and the obtained material is used to assess compound activity in an autophosphorylation assay. Purified enzyme (10 nM) is incubated in 50 μL of Tris buffer (Tris 50 mM, pH 7.4; NaCl, 100 mM; MgCl ₂ , 5 mM; MnCl ₂ , 5 mM; and DTT, 10 mM). The enzyme is preincubated with different concentrations of RepSox (0.1% DMSO final concentration in the test) for 10 min at 37°C. The reaction is then initiated by the addition of 3 μM ATP (0.5 μCi γ-33P-ATP). After 15 min at 37°C, phosphorylation is stopped by the addition of SDS-PAGE sample buffer (50 mM Tris-HCl,

A DRUG SCREENING EXPERT

Kinase Assay	pH 6.9, 2.5% glycerol, 1% SDS, and 5% β -mercaptoethanol). The samples are boiled for 5 min at 95°C and run on a 12% SDS-PAGE. Dried gels are exposed to a phosphor screen overnight. ALK5 autophosphorylation is quantified using a Storm imaging system[1].
Cell Research	RepSox is dissolved in DMSO and stored, and then diluted with appropriate media (DMSO 1%) before use[1]. To test anti-TGF- β activity of compounds, HepG2 cells are seeded in 96 well microplates at a concentration of 35000 cells per well in 200 μ L of serum-containing medium. The microplates are then placed for 24 h in a cell incubator at 37°C, 5% CO ₂ atm. RepSox dissolved in DMSO are then added at concentrations of 50 nM to 10 μ M (final concentration of DMSO 1%) for 30 min prior to the addition of recombinant TGF- β (1 ng/mL). After an overnight incubation, the cells are washed with PBS and lysed by addition of 10 μ L of passive lysis buffer. Inhibition of luciferase activity relative to control groups is used as a measure of compound activity. A concentration-response curve is constructed from which an IC ₅₀ value is determined graphically[1].

Solubility Information

Solubility	Ethanol: < 1 mg/mL (insoluble or slightly soluble), H ₂ O: < 1 mg/mL (insoluble or slightly soluble), DMSO: 51 mg/mL (177.5 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 5 mg/mL (17.4 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4804 mL	17.4022 mL	34.8044 mL
5 mM	0.6961 mL	3.4804 mL	6.9609 mL
10 mM	0.348 mL	1.7402 mL	3.4804 mL
50 mM	0.0696 mL	0.348 mL	0.6961 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Tu WZ, et al. RepSox, a small molecule inhibitor of the TGF β receptor, induces brown adipogenesis and browning of white adipocytes. *Acta Pharmacol Sin.* 2019 Dec;40(12):1523-1531.
- Ma X, Lu Y, Zhou Z, Human expandable pancreatic progenitor-derived β cells ameliorate diabetes. *Science Advances.* 2022, 8(8): eabk1826.
- Wang W, Ren S, Lu Y, et al. Inhibition of Syk promotes chemical reprogramming of fibroblasts via metabolic rewiring and H2S production. *The EMBO Journal.* 2021 Jun 1;40(11):e106771. doi: 10.15252/embj.2020106771. Epub 2021 Apr 28.
- Guo Y, et al. RepSox effectively promotes the induced differentiation of sheep fibroblasts into adipocytes via the inhibition of the TGF- β 1/Smad pathway. *Int J Mol Med.* 2021 Aug;48(2):148.
- Shi CJ, et al. TGF β R-1/ALK5 inhibitor RepSox induces enteric glia-to-neuron transition and influences gastrointestinal mobility in adult mice. *Acta Pharmacol Sin.* 2023 Jan;44(1):92-104.
- Li J, Bai Y, Liu Y, et al. Transcriptome-based chemical screens identify CDK8 as a common barrier in multiple cell reprogramming systems. *Cell Reports.* 2023, 42(6).
- Zhang S, Lou H, Lu H, et al. Characterization of Proliferation Medium and Its Effect on Differentiation of Muscle Satellite Cells from *Larimichthys crocea* in Cultured Fish Meat Production. *Fishes.* 2023, 8(9): 429.
- He D, et al. TGF- β inhibitor RepSox suppresses osteosarcoma via the JNK/Smad3 signaling pathway. *Int J Oncol.* 2021 Nov;59(5):84.
- Xin Q, Niu R, Chen Q, et al. Stable cytoactivity of piscine satellite cells in rice bran-gelatin hydrogel scaffold of cultured meat. *International Journal of Biological Macromolecules.* 2024: 134242.
- Lou H, Lu H, Zhang S, et al. Highly aligned myotubes formation of piscine satellite cells in 3D fibrin hydrogels of cultured meat. *International Journal of Biological Macromolecules.* 2024: 136879.
- Zhang S, Lu H, Lou H, et al. An efficient serum-free medium for ex vivo expansion of myoblasts from *Larimichthys crocea* for cultured meat production. *Food Research International.* 2024, 196: 115073.
- Qin X, Fu L, Li C, et al. Optimized inner ear organoids for efficient hair cell generation and ototoxicity response modeling. *Science China Life Sciences.* 2025: 1-15.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481