

Anti-HER3/ERBB3 Antibody (2Y512)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-HER3/ERBB3 Antibody (2Y512) is a mouse (variable region) / human (kappa / IgG1 constant) chimeric monoclonal antibody targeting Human HER3/ERBB3.
Targets(IC50)	HER

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481