

Anti-TNFR1/CD120a/TNFRSF1A Antibody (9Z795)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|------------------------------------------------------------------------------------------------------------------------|
| Description | Anti-TNFR1/CD120a/TNFRSF1A Antibody (9Z795) is a Rabbit IgG monoclonal antibody targeting Human TNFR1/CD120a/TNFRSF1A. |
| Targets(IC50) | TNF |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481