

Anti-Human respiratory syncytial virus (RSV) Fusion Glycoprotein/RSV-F Antibody (1B191)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Human respiratory syncytial virus (RSV) Fusion Glycoprotein/RSV-F Antibody (1B191) is a Rabbit IgG monoclonal antibody targeting RSV Human respiratory syncytial virus (RSV) Fusion Glycoprotein/RSV-F.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481