

## Anti-VEGF165/VEGFA Antibody (5L924)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-VEGF165/VEGFA Antibody (5L924) is a Rabbit IgG monoclonal antibody targeting Human, Cynomolgus VEGF165/VEGFA.
Targets(IC50)	VEGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481